



INCYT RESILIENCE INC.

Sold By Tri Hawk Canada
www.trihawk.com



Made in Canada

We manufacture in Canada. We test in Canada.

We are a team of Canadians that are building new capabilities and adding the capacity for manufacturing that Canada needs - today, and in the future. We believe Canada can build and manufacture with the best in the world.

Shop online now!

<https://www.trihawk.com>



Contact us

Tri Hawk Canada/Xenopus Inc.
1-877—TRI HAWK (874-4295)
+1-514-866-4295
dentalinfo@trihawk.com

ASTM Level 3 Surgical Mask

Medical mask test summary
Tested in Canada by Kinectrics



Particulate Filtration Efficiency
ASTM F2299

> 98%
(99.93 %)

Bacterial Filtration Efficiency
ASTM F 2101

> 98%
(99.76 %)

Differential Pressure
EN 14683 Annex C

< 6.0 mm H₂O/cm²
(4.41 mm H₂O/cm²)

Fluid Resistance to Synthetic Blood
ASTM F1862

> 160 mmHg
(pass)

Flammability
16 CFR 1610

Class 1



* Available in 50-mask and 10-mask retail boxes. Cases of 25 boxes also available.